


Product Datasheet

Murine FGF-8 protein (orb107713)

Description	Fibroblast Growth Factor-8 (FGF-8) is a heparin binding growth factor belonging to the FGF family. P...
Conjugation	Unconjugated
Tested Applications	HPLC, SDS-PAGE
Target	FGF-8
Preservatives	Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH 7.4 500mMNaCl.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	QVRSAAQKRG PGAGNPADTL GQGHEDRPFG QRSRAGKNFT NPAPNYPEEGSKEQRDSVLP KVTQRHVREQ SLVTDQLSRR LIRTYQLYSR TSGKHVQVLA NKRINAMAEDGDPFAKLIVE TDTFGSRVRV RGAETGLYIC MNKKGKLIK SNGKGKDCVF TEIVLENNYTALQNAKYEGW YMAFTRKGRP RKGSKTRQHQ REVHFMKRLP RGHHTTEQSLRFEFLNYPPF TRSLRGSQRT WAPEPR
Purity	> 97% by SDS-PAGE and HPLC analyses
MW	28.1kDa
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
Murine FGF-8
protei...