


Product Datasheet

Rat FGF-10 protein (orb107717)

Description	FGF-10 was originally identified from rat embryos by homology-based polymerase chain reaction. Rat F...
Conjugation	Unconjugated
Tested Applications	HPLC, SDS-PAGE
Target	FGF-10
Preservatives	Lyophilized from a 0.2 µm filtered concentrated solution in 20mM Tris, 500mM NaCl, pH7.4
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	QALGQDMVSP EATNSSSSSS SSSSSSFSS PSSAGRHVRS YNHLQGDVRW RKLFSFTKYFLKIEKNGKVS GTKKENCYPYS ILEITSVEIG VVAVKAINSN YYLAMNKKGK LYGSKEFNNDCKLKERIEEN GYNTYASFNW QHNGRQMYVA LNGKGAPRRG QKTRRKNTSAHFLPMVVHS
Purity	> 96% by SDS-PAGE and HPLC analyses
MW	20.0 kDa
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
Rat FGF-10
protein