


Product Datasheet

Murine FGF-21 protein (orb107722)

Description	Recombinant Murine Fibroblast Growth Factor-21. Fibroblast growth factor-21 (FGF-21) belongs to the l...
Conjugation	Unconjugated
Tested Applications	HPLC, SDS-PAGE
Target	FGF-21
Preservatives	Lyophilized from a 0.2 µm filtered concentrated solution in 3xPBS, pH 7.4
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	AYPIPDSSPL LQFGGQVRQR YLYTDDDQDT EAHLEIREDG TVVGAHRSP ESLLLELKALPGVIQILGVK ASRFLCQQPD GALYGSPHFD PEACSFRELL LEDGYNVYQS EAHGLPLRLPQKDSPNQDAT SWGPVRFLPM PGLLHEPQDQ AGFLPPEPPD VGSSDPLSMV EPLQGRSPSYAS
Purity	> 97% by SDS-PAGE and HPLC analyses
MW	19.9 kDa
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Murine FGF-21 prote...