

Product Datasheet

Human FGF-23 protein (orb107724)

Description	Fibroblast growth factor-23 (FGF-23) belongs to the large FGF family which has at least 23 members. ...
Conjugation	Unconjugated
Tested Applications	HPLC, SDS-PAGE
Target	FGF-23
Preservatives	Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH 7.4.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	MYPNASPLLGSWGGGLIHLYTATARNSYHL QIHKNGHVDPGAPHQTIYSALMIRSEDAGFV VITGVMSRRYL CMDFRGNIFGSHYFDPENC RFQHQTLENGYDVYHSPQYHFLVSLGRAKR AFLPGMNPPPYSQFLSRRNE IPLIHFNTPIPRRHTRSAEDDSERDPLNLV KPRARMTPAPASCSQELPSA EDNSPMASDPLGVVRGGRVNTHAGGTGPEG CRPFAKFI
Purity	> 95% by SDS-PAGE and HPLC analyses
MW	22.5 kDa
Expiration Date	6 months from date of receipt.