


Product Datasheet

Rat Fractalkine protein (orb107726)

Description	Recombinant Rat Fractalkine/CX3CL1. CX3CL1 recently identified through bioinformatics is the only kn...
Conjugation	Unconjugated
Tested Applications	HPLC, SDS-PAGE
Target	Fractalkine
Preservatives	Lyophilized from a 0.2 µm filtered concentrated solution in 2xPBS, pH 7.4.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	QHLGMTKCNI TCHKMTSPIP VTLLIHYNLN QESCGKRAII LETRQHRHFC ADPKEKWVDAMKHLDHQTA ALTRNG
Purity	> 95% by SDS-PAGE and HPLC analyses
MW	8.8 kDa
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
Rat
Fractalkine
pro...