


Product Datasheet

Human GDNF protein (orb107727)

Description	Glial cell-derived neurotrophic factor, also known as GDNF is a small protein that potently promotes...
Conjugation	Unconjugated
Tested Applications	HPLC, SDS-PAGE
Target	GDNF
Preservatives	Lyophilized from a 0.2 µm filtered concentrated solution in 20mM PB, 500mM NaCl, pH 6.4
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	MSPDKQMAVL PRRERNRQAA AANPENSRGK GRRGQRGKNR GCVLTAIHLNVTDLGLGYET KEELIFRYCS GSCDAAETTY DKILKNLSRN RRLVSDKVGQ ACCRPIAFDDDL SFLDDNLV YHILRKHS AK RCGCI
Purity	> 97% by SDS-PAGE and HPLC analyses.
MW	15.2 kDa
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human GDNF protein