


Product Datasheet

Human HB-EGF protein (orb107732)

Description	Human HB-EGF protein
Conjugation	Unconjugated
Tested Applications	HPLC, SDS-PAGE
Target	Human HB-EGF
Preservatives	Lyophilized from a 0.2 µm filtered concentrated solution in 20mM PB, pH7.4, 130mM NaCl
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	DLQEADLDLL RVTLSKPKQA LATPNKEEHG KRKKKGKGLG KKRDPCLRKYKDFCIHGECK YVKELRAPSC ICHPGYHGER CHGLSL
Purity	> 97% by SDS-PAGE and HPLC analyses
MW	9.7 kDa
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human HB-EGF protei...