

Product Datasheet

HBsAg antibody (orb10774)

Description	HBsAg antibody
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant HBsAg protein
Target	HBsAg
Preservatives	0.01M PBS(pH7.4) with 0.03% Proclin300
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgM
Clonality	Monoclonal
MW	44 kDa
Uniprot ID	Q89714
Expiration Date	12 months from date of receipt.