


Product Datasheet

Rat IL1 alpha protein (orb107740)

Description	Recombinant rat IL1 alpha protein
Conjugation	Unconjugated
Tested Applications	HPLC, SDS-PAGE
Target	IL1 alpha
Preservatives	Lyophilized from a 0.2 µm filtered concentrated solution in 20 mM Tris, 150 mM NaCl, pH 8.5
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	SAPHSFQNNL RYKLIRIVKQ EFIMNDSLQ NIYVDMRIH LKAASLNDLQLEVKFDYAY SSGGDDSKYP VTLKVSNTQL FVSAQGEDKP VLLKEIPETPKLITGSETDL IFFWEKINSK NYFTSAAFPE LLIATKEQSQ VHLARGLPMSIDFQIS
Purity	> 97 % by SDS-PAGE and HPLC analyses
MW	17.8 kDa
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
Rat IL1 alpha
prote...