


Product Datasheet

Rat IL4 protein (orb107746)

Description	Rat IL4 protein
Conjugation	Unconjugated
Tested Applications	HPLC, SDS-PAGE
Target	Rat IL4
Preservatives	Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH 7.4, 5% trehalose.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	HGCNDSPLRE IINTLNQVTE KGTPCTEMFV PDVLTATRNT TENELICRAS RVLRFYFPRDVPPCLKNKS GVLGELRKLC RGVSGLSLR SCTVNESTLT TLKDFLES�K SILRGKYLQSCTSMS
Purity	> 95 % by SDS-PAGE and HPLC analyses
MW	14.0 kDa
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
Rat IL4
protein