


Product Datasheet

Rhesus IL4 protein (orb107747)

Description	Rhesus IL4 protein
Conjugation	Unconjugated
Tested Applications	HPLC, SDS-PAGE
Target	Rhesus IL4
Preservatives	Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH 7.4, 3 % trehalose
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	HNCHIALREI IETLNSLTEQ KTLCTKLTIT DILAASKNTT EKETFCRAAT VLRQFYSHHEKDTRCLGATA QQFHRHKQLI RFLKRLDRNL WGLAGLNSCP VKEANQSTLE DFLERLKTIMREKYSKCSS
Purity	> 96 % by SDS-PAGE and HPLC analyses
MW	15.0 kDa
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Rhesus IL4 protein