


Product Datasheet

Human IL-7 protein (orb107750)

Description	Interleukin-7 (IL-7) is encoded by the IL7 gene and secreted by stromal cells in the red marrow and ...
Conjugation	Unconjugated
Tested Applications	HPLC, SDS-PAGE
Target	IL-7
Preservatives	Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH 7.4, 2% trehalose
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	DCDIEGKDG KQYESVLMVS IDQLLDSTMKE IGSNCLNNEF NFFKRHICDANKEGMFLFRA ARKLRQFLKM NSTGDFDLHL LKVSEGTIL LNCTGQVKGRKPAALGEAQP TKSLEENKSL KEQKKLNDLC FLKRLLEIK TCWNKILMGT KEH
Purity	> 96% by SDS-PAGE and HPLC analyses
MW	17.4 kDa
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human IL-7 protein