


Product Datasheet

Rat IL-7 protein (orb107751)

Description	Interleukin-7 (IL-7) is encoded by the IL7 gene and secreted by stromal cells in the red marrow and ...
Conjugation	Unconjugated
Tested Applications	HPLC, SDS-PAGE
Target	IL-7
Preservatives	Lyophilized from a 0.2 m filtered concentrated solution in PBS, pH 7.4.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	DCHIKDKDGK AFGSVLMISI NQLDKMTGTD SDCPNNEPNF FKKHLCDDTKEAAFLNRAAR KLRQFLKMNI SEEFNDHLLR VSDGTQTLVN CTSKEEKTIKEQKNDPCFL KRLLEIKTC WNKILKGS
Purity	> 95 % by SDS-PAGE and HPLC analyses
MW	14.9 kDa
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
Rat IL-7
protein