


Product Datasheet

Rhesus IL-8 protein (orb107752)

Description	Interleukin-8 (C-X-C motif chemokine 8) is coded by the interleukin 8 gene on the chromosome 5 in Rh...
Conjugation	Unconjugated
Tested Applications	HPLC, SDS-PAGE
Target	IL-8
Preservatives	Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH 7.4
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	AVLPRSAKEL RCECIKTYSK PFHPKFIKEL RVIESGPHCA NTEIIVKLSD GRELCCLDPKEPWVQRVVEKF VKRAENQNP
Purity	> 95% by SDS-PAGE and HPLC analyses
MW	9.14kD
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Rhesus IL-8 protein