


Product Datasheet

Murine IL-16 protein (orb107757)

Description	Interleukin-16 (IL-16) is also called lymphocyte chemoattractant factor (LCF) and it is mostly secre...
Conjugation	Unconjugated
Tested Applications	HPLC, SDS-PAGE
Target	IL-16
Preservatives	Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH 7.4
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	SAASASAASD ISVESKEATV CTVTLEKTS A GLGFSLEGGK GSLHGDKPLT INRIFKGTEQGEMVQPGDEI LQLAGTAVQG LTRFEAWNVI KALPDGPVTI VIRRTSLQCK QTTASADS
Purity	> 98 % by SDS-PAGE and HPLC analyses.
MW	12.2 kDa
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Murine IL-16 protei...