

Product Datasheet

Rat IL-33 protein (orb107761)

Description	Interleukin-33 (IL-33), also known as NF-HEV and DVS 27, is a cytokine belonging to the IL-1 superfa...
Conjugation	Unconjugated
Tested Applications	HPLC, SDS-PAGE
Target	IL-33
Preservatives	Lyophilized from a 0.2 µm filtered concentrated solution in 20 mM Tris, 300 mM NaCl, pH 8.5, 3mM beta -mercaptoethanol.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	SIQGTSLLTE SCALSTYNDQ SVSFVLENGC YVINVEDCGK NQEKDKVLLR YYESSFPAQSGDGVGKKLM VNMSPIKDTD IWLNANDKDY SVELQKGDVS PPDQAFFVLHKKSSDFVSFE CKNLPGTYIG VKDNQLALVE ENDESCNNIM FKLSKM
Purity	> 95 % by SDS-PAGE and HPLC analyses
MW	17.4 kDa
Expiration Date	6 months from date of receipt.