


Product Datasheet

Human IL-36RA protein (orb107770)

Description	Interleukin-36 receptor antagonist (IL-36RA) is a secreted protein which belongs to the interleukin ...
Conjugation	Unconjugated
Tested Applications	HPLC, SDS-PAGE
Target	IL-36RA
Preservatives	Lyophilized from a 0.2 "m filtered concentrated solution in PBS, pH 7.4.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	VLSGALCFRM KDSALKVLYL HNNQLLAGGL HAGKVIKGE ISVVPNRWLD ASLSPVILGVQGSQCLSCG VGQEPTLTLE PVNIMELYLG AKESKSFTFY RRDMLTSSF ESAAYPGWFLCTVPEADQPV RLTQLPENGW WNAPITDFYF QQCD
Purity	> 95 % by SDS-PAGE and HPLC analyses
MW	16.8 kDa
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
Human IL-
36RA prote...