

Product Datasheet

Murine IL-36RA protein (orb107771)

Description	Interleukin-36 receptor antagonist (IL-36RA) is a secreted protein which belongs to the interleukin ...
Conjugation	Unconjugated
Tested Applications	HPLC, SDS-PAGE
Target	IL-36RA
Preservatives	Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH 7.4, 2mM beta -mercaptoethanol
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	VLSGALCFRM KDSALKVLYL HNNQLLAGGL HAEKVIKGEE ISVVPNRALD ASLSPVILGVQGGSQLSCG TEKGPILKLE PVNIMELYLG AKESKSFTFY RRDMGLTSSF ESAAYPGWFLCTSPEADQPV RLTQIPEDPA WDAPITDFYF QQCD
Purity	> 95 % by SDS-PAGE and HPLC analyses
MW	16.9 kDa
Expiration Date	6 months from date of receipt.