

Product Datasheet

Murine IP10 protein (orb107772)

Description	Recombinant Murine gamma-Interferon Inducible Protein 10/CXCL10. CXCL10 also known as IP-10 is belon...
Conjugation	Unconjugated
Tested Applications	HPLC, SDS-PAGE
Target	IP10
Preservatives	Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH 7.4.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	IPLARTVRCN CIHIDDPVR MRAIGKLEII PASLSCPVE IIATMKKNDE QRCLNPESKTIKMLKAFSQ KRSKRAP
Purity	> 97% by SDS-PAGE and HPLC analyses.
MW	8.7 kDa
Expiration Date	6 months from date of receipt.