


Product Datasheet

Rat IP-10 protein (orb107773)

Description	gamma -Interferon Inducible Protein 10 (IP-10)/CXCL10 was originally identified as an IFN-gamma-indu...
Conjugation	Unconjugated
Tested Applications	HPLC, SDS-PAGE
Target	IP-10
Preservatives	Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH 7.4.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	IPLARTVRCT CIDFHEQPLR PRAIGKLEII PASLSCPHVE IIATMKKNNE KRCLNPESEAIKSLKAVSQ RRSKRAP
Purity	> 95% by SDS-PAGE and HPLC analyses
MW	8.6 kDa
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Rat IP-10 protein