


Product Datasheet

Rat LIX protein (orb107775)

Description	Recombinant Rat LPS-induced CXC chemokine/CXCL5. CXCL5 is a small cytokine belonging to the CXC chem...
Conjugation	Unconjugated
Tested Applications	HPLC, SDS-PAGE
Target	LIX
Preservatives	Lyophilized from a 0.2 µm filtered concentrated solution in 2xPBS, pH 7.4.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	APSSVIAATE LRCVCLTVTP KINPKLIANL EVIPAGPQCP TVEVIAKLKNQKEVCLDPEA PVIKKIIIQK ILGSDKKKAK RNALAVERTA SVQ
Purity	> 97% by SDS-PAGE and HPLC analyses
MW	9.8 kDa
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Rat LIX protein