


Product Datasheet

Human MANF protein (orb107776)

Description	Mesencephalic astrocyte-derived neurotrophic factor, also known as Protein ARMET, Arginine-rich prot...
Conjugation	Unconjugated
Tested Applications	HPLC, SDS-PAGE
Target	MANF
Preservatives	Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH 7.4.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	LRPGDCEVCI SYLGRFYQDL KDRDVTFSPA TIENELIKFC REARGKENRL CYYIGATDDAATKIINEVSK PLAHHIPVEK ICEKLKKKDS QICELKYDKQ IDLSTVDLKK LRVKELKKILDDWGETCKGC AEKSDYIRKI NELMPKYAPK AASARTDL
Purity	> 95% by SDS-PAGE and HPLC analyses.
MW	18.2 kDa
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
Human MANF
protein