

Product Datasheet

Murine M-CSF protein (orb107778)

Description	Macrophage Colony Stimulating Factor (M-CSF), also known as CSF-1, is a four F helical bundle cytoki...
Conjugation	Unconjugated
Tested Applications	HPLC, SDS-PAGE
Target	M-CSF
Preservatives	Lyophilized from a 0.2 µm filtered solution in 20mM Tris, 500mM NaCl, pH7.4.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	KEVSEHCSTM IGNGHLKVLQ QLIDSQMETS CQIAFEFVDQ EQLDDPVCYL KKAFFLVQDIIDETMRFKDN TPANATERL QELSNNLNSC FTKDYEEQNK ACVRTFHETP LQLEKIKNFFNETKNLLEK DWNIFTKNCN NSFAKCSSRD VVTKPDCNCL YPKATPSSDP ASASPHQPPAPSMAPLAGLA WDDSQRTEGS SLLPSELPLR IEDPGSAKQR PPRSTCQTLE
Purity	> 95% by SDS-PAGE and HPLC analyses.
MW	25.9kDa
Expiration Date	6 months from date of receipt.