


Product Datasheet

Murine MCP3 protein (orb107779)

Description	Murine MCP3 protein
Conjugation	Unconjugated
Tested Applications	HPLC, SDS-PAGE
Target	Murine MCP3
Preservatives	Lyophilized from a 0.2 µm filtered concentrated solution in 2xPBS, pH 7.4.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	QPDGPNASTC CYVKKQKIPK RNLKSYRRIT SSRCPWEAVI FKTKKGMEVCAEAHQKWVEE AIAYLDMKTP TPKP
Purity	> 95% by SDS-PAGE and HPLC analyses
MW	8.5 kDa
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
Murine MCP3
protein