


Product Datasheet

Rat MEC protein (orb107783)

Description	Mucosae-associated Epithelial Chemokine (MEC)/CCL28 (CC chemokine ligand 28) is a secreted CC chemok...
Conjugation	Unconjugated
Tested Applications	HPLC, SDS-PAGE
Target	MEC
Preservatives	Lyophilized from a 0.2 µm filtered concentrated solution in 20mM PB, pH 7.4, 200mM NaCl.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	SEAILPIASS CTEVSHHIP RLLERVNSC SIQRADGDCD LAAVILHVKR RRICVSPHNPTLKRWMSASE MKNGKENLCP RKKQDSGKDR KGHTPRKHGK HGTRRIHGTH DHEAPR
Purity	> 96% by SDS-PAGE and HPLC analyses.
MW	13.1 kDa
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Rat MEC protein