

Product Datasheet

Rat Oncostatin M protein (orb107792)

Description	Recombinant Rat Oncostatin-M (OSM) is a multifunctional cytokine that belongs to the Interleukin-6 s...
Conjugation	Unconjugated
Tested Applications	HPLC, SDS-PAGE
Target	Oncostatin M
Preservatives	Lyophilized from a 0.2 µm filtered concentrated solution in 2xPBS, 1 % Tween-80
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	KRGCSSSSPK LLSQLKSQAN ITGNTASLLE PYILHQNLNT LTLRAACTEH PVAFPSEDMLRQLSKPDFLS TVHATLGRVW HQLGAFRQQF PKIQDFPELE RARQNIQGIR NNVYCMARLLHPPLEIPEPT QADSGTSRPT TTAPGIFQIK IDSCRFLWGY HRFMGSVGRV FEEWGDGSRRSRRHSPLWAW LKGDHRIRPS RSSQSAMLRS LVPR
Purity	> 96 % by SDS-PAGE and HPLC analyses
MW	24.3 kDa
Expiration Date	6 months from date of receipt.