

Product Datasheet

Rhesus S100B protein (orb107799)

Description	Rhesus S100B protein
Conjugation	Unconjugated
Tested Applications	HPLC, SDS-PAGE
Target	S100B
Preservatives	Lyophilized from a 0.2 µm filtered solution in PBS, pH 7.4
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	SELEKAMVAL IDVFHQYSGR EGDKHKLKKS ELKELINNEL SHFLEEIQ EVVDKVMETLSDGDGECDF QEFMAFVAMV TTACHEFFE E
Purity	> 97% by SDS-PAGE and HPLC analyses
MW	10.7 kDa
Expiration Date	6 months from date of receipt.