

Product Datasheet

Rat SCF protein (orb107801)

Description	Stem Cell Factor (SCF) which binds to the c-Kit receptor is produced by fibroblasts and endothelial ...
Conjugation	Unconjugated
Tested Applications	HPLC, SDS-PAGE
Target	SCF
Preservatives	Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH 7.4
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	MQEICRNPVT DNVKDITKLV ANLPNDYMIT LNYVAGMDVL PSHCWLRDMVTHLSVSLTTL LDKFSNISEG LSNYSIIDKL GKIVDDLVA MEENAPKNVK ESLKKPETRNFTPEEFSIF NRSIDAFKDF MVASDTSDCV LSSTLGPEKD SRVSVTKPFM LPPVA
Purity	> 95% by SDS-PAGE and HPLC analyses
MW	18.5 kDa
Expiration Date	6 months from date of receipt.