


Product Datasheet

Murine SHH protein (orb107802)

Description	Sonic hedgehog is one of three proteins in the mammalian signaling pathway family called hedgehog, t...
Conjugation	Unconjugated
Tested Applications	HPLC, SDS-PAGE
Target	SHH
Preservatives	Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH 7.4.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	IVIGPGRGFG KRRHPKKLTP LAYKQFIPNV AEKTLGASGR YEGKITRNSERFKELTPNYN PDIIFKDEEN TGADRLMTQR CKDKLNALAI SVMNQWPGVKLRVTEGWDED GHHSEESLHY EGRAVDITTS DRDRSKYGML ARLAVEAGFDWVYYESKAHI HCSVKAENSV AAKSGG
Purity	> 95 % by SDS-PAGE and HPLC analyses
MW	19.8 kDa
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
Murine SHH
protein