

Product Datasheet

Human SYND4 protein (orb107803)

Description	Syndecan-4 Human Recombinant fused with 6X His tag produced in E.coli is a single, non-glycosylated,...
Conjugation	Unconjugated
Tested Applications	HPLC, SDS-PAGE
Target	SYND4
Preservatives	Lyophilized from a 0.2 µm filtered concentrated solution in 20mM PB, pH 7.4.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Purity	> 90% by SDS-PAGE and HPLC analyses
MW	15.4 kDa
Uniprot ID	P31431
Solubility (25°C)	It is recommended to reconstitute the lyophilized SDC4 in sterile H ₂ O not less than 100 ug/ml, which can then be further diluted to other aqueous solutions.
NCBI	NM_002999.3, NP_002990.2
Entrez	6385
Expiration Date	6 months from date of receipt.