


Product Datasheet

Murine TECK protein (orb107807)

Description	Thymus Expressed Chemokine (TECK)/CCL25 is a novel CC chemokine that is distantly related (approxima...
Conjugation	Unconjugated
Tested Applications	HPLC, SDS-PAGE
Target	TECK
Preservatives	Lyophilized from a 0.2 µm filtered concentrated solution in 20mM PB, pH 7.4, 150mM NaCl
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	QGAFEDCCLG YQHRIKWNVL RHARNYHQQE VSGSCNLRAV RFYFRQKVVCGNPEDMNVKR AIRILTARKR LVHWKSASDS QTERKKSNNH KSKVENPNSTSVRSATLGHP RMVMMPRKTN N
Purity	> 95% by SDS-PAGE and HPLC analyses.
MW	14.1 kDa
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
Murine TECK
protein