


Product Datasheet

Human sTRAIL Receptor-2 protein (orb107809)

Description	Tumor necrosis factor-related apoptosis-inducing ligand Receptor 2 (TRAIL-R2) is a cell-surface rece...
Conjugation	Unconjugated
Tested Applications	HPLC, SDS-PAGE
Target	sTRAIL Receptor-2
Preservatives	Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH 7.4.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	ESALITQQDL APQQR AAPQ KRSSPSEGLC PPGHHISEDG RDCISCKYGQ DYSTHWNDLLFCLRCTRCDS GEVELSPCTT TRNTVCQCEE GTFREEDSPE MCRKCRTGCP RGMVKVGDC TPWSDIECVHK ES
Purity	> 97% by SDS-PAGE and HPLC analyses
MW	14.8 kDa
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human sTRAIL Recept...