


Product Datasheet

Human UBC1 protein (orb107811)

Description	Ubiquitin-conjugating enzyme E2 K is a protein that in humans is encoded by the UBE2K gene. The prot...
Conjugation	Unconjugated
Tested Applications	HPLC, SDS-PAGE
Target	UBC1
Preservatives	Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH 7.4.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	MHHHHHHHAMA NIAVQRIKRE FKEVLKSEET SKNQIKVDLV DENFTELRGEIAGPPDTPYE GGRYQLEIKI PETYPFNPPK VRFITKIWHP NISSVTGAICLDILKDQWAA AMTLRTVLLS LQALLAAEP DDPQDAVVAN QYKQNPENFKQTARLWAHVY AGAPVSSPEY TKKIENLCAM GFDRNAVIVA LSSKSWDVETATELLLSN
Purity	> 95% by SDS-PAGE and HPLC analyses
MW	23.4 kDa
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
Human UBC1
protein