


## Product Datasheet

### Human UBC2 protein (orb107812)

<b>Description</b>	This E2 enzyme encodes for the human homolog of the yeast DNA repair gene RAD6, which is induced by ...
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	HPLC, SDS-PAGE
<b>Target</b>	UBC2
<b>Preservatives</b>	Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH 7.4
<b>Storage</b>	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Protein Sequence</b>	MHHHHHHAMG QLRSMSTPAR RRLMRDFKRL QEDPPVGVS APSENNIMQWNAVIFGPEGT PFEDGTFKLV IEFSEEYPNK PPTVRFLSKM FHPNVYADGS ICLDILQNRWSPTYDVSSIL TSIQSLLDEP NPNSPANSQA AQLYQENKRE YEKRVSVAIVE QSWNDS
<b>Purity</b>	> 95% by SDS-PAGE and HPLC analyses
<b>MW</b>	19.0 kDa
<b>Expiration Date</b>	6 months from date of receipt.



SDS-PAGE analysis of Human UBC2 protein