

Product Datasheet

Human UBC7 protein (orb107814)

Description	Human Ubiquitin Conjugating Enzyme E2 L3 (UBE2L3)/UBCH7 is a class I enzyme which functions in the st...
Conjugation	Unconjugated
Tested Applications	HPLC, SDS-PAGE
Target	UBC7
Preservatives	Lyophilized from a 0.2 µm filtered concentrated solution in PBS, 1mM DTT, pH 7.5
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	MHHHHHHHAMA ASRRLMKELE EIRKCGMKNF RNIQVDEANL LTWQGLIVPDNPPYDKGAFR IEINFPAEYP FKPPKITFKT KIYHPNIDEK GQVCLPVISA ENWKPATKTDQVIQSLIALV NDPQPEHPLR ADLAE EYSKD RKKFCKNAEE FTKKYGEKRP VD
Purity	> 95% by SDS-PAGE and HPLC analyses.
MW	18.9 kDa,
Expiration Date	6 months from date of receipt.