

Product Datasheet

Rat VEGF165 protein (orb107815)

Description	Vascular Endothelial Growth Factor is a potent growth and angiogenic cytokine. It stimulates prolif...
Conjugation	Unconjugated
Tested Applications	HPLC, SDS-PAGE
Target	VEGF165
Preservatives	Lyophilized from a 0.2 µm filtered concentrated solution in 20 mM Tris, 300 mM NaCl, pH 8.8.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	MAPTTEGEQK AHEVVKFMDV YQRSYCRPIE TLVDIFQEYP DEIEYIFKPS CVPLMRCAGCCNDEALECVP TSESNVTMQI MRIKPHQSQH IGEMSFLQHS RCECRPKKDR TKPEKHCEPCSERRKHLFVQ DPQTCKCSCK NTDSRCKARQ LELNERTCRC DKPRR
Purity	> 95 % by SDS-PAGE and HPLC analyses
MW	19.4 kDa
Expiration Date	6 months from date of receipt.