



Product Datasheet

HDAC5 antibody (orb10796)


Description	Rabbit polyclonal antibody to HDAC5
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	A synthetic peptide of human HDAC5(NP_005465.2).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	132kDa
Uniprot ID	Q9UQL6
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of A-549 cells.



Immunohistochemical staining of mouse spleen.



Immunohistochemical staining of human tissue.