

Product Datasheet

HHV8 antibody (orb10800)

Description	HHV8 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	KLH conjugated synthetic peptide derived from human HHV8/ORF50 (501-600/691 aa)
Target	HHV8 ORF50
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	74 kDa
Uniprot ID	F5HCV3
Expiration Date	12 months from date of receipt.