

Product Datasheet

Anti-beta-Actin/Cy5.5 (orb1081762)

Description	Anti-beta-Actin/Cy5.5
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Cy5.5
Tested Applications	FC, ICC, IF
Immunogen	Synthetic MAP peptide derived from human beta-Actin
Target	beta-Actin (Loading Control)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P60709
Dilution Range	Flow-Cyt=1ug/Test,IF=1:200-500,ICC=1:100
Expiration Date	12 months from date of receipt.