

Product Datasheet

HIV gp41 antibody (orb10820)

Description	HIV gp41 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	KLH conjugated synthetic peptide derived from Transmembrane protein gp41 (1-80/120 aa)
Target	HIV gp41
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	39/93 kDa
Uniprot ID	Q79670
Expiration Date	12 months from date of receipt.