



Product Datasheet

HLA-A antibody (orb10824)

Description	HLA-A antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from the middle of human HLA-G (181-280/338 aa)
Target	HLA G
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	35 kDa
Uniprot ID	P17693
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



IHC-P of human lung cancer tissue (HLA G...)



Western blot analysis of mouse brain lys...