


Product Datasheet

HNF 1 alpha antibody (orb10829)


Description	HNF 1 alpha antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human HNF1 (201-350/628 aa)
Target	HNF1A
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	67 kDa
Uniprot ID	P15257
Dilution Range	WB=1:500-2000, IHC-P=1:100-500, ICC=1:100, Flow-Cyt=3ug/test
Expiration Date	12 months from date of receipt.



Western blot analysis
of Human HepG2
Cel...



Immunohistochemical
staining of Mouse li...



Immunohistochemical
staining of Rat live...