

Product Datasheet

Anti-phospho-C-Myc(Ser62)/Cy5.5 [Out of stock] (orb1082971)

Description	Rabbit polyclonal antibody to Anti-phospho-C-Myc(Ser62)/Cy5.5
Species/Host	Rabbit
Conjugation	Cy5.5
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human C-Myc around the phosphorylation site of Ser62
Preservatives	Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4.
Form/Appearance	Lyophilized or Liquid
Concentration	1mg/ml
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	49kDa
Uniprot ID	P01106, P01108, P09416
Expiration Date	12 months from date of receipt.