

Product Datasheet

HP1 alpha (phospho-Tyr24) antibody (orb10834)

Description	HP1 alpha (phospho-Tyr24) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from Fruit Fly HP1 around the phosphorylation site of Tyr24 (EE(p-Y)AV)
Target	phospho-heterochromatin protein 1 (pTyr24)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	23 kDa
Uniprot ID	P05205
Expiration Date	12 months from date of receipt.