



Product Datasheet

HPV18-E6/E6 antibody (orb10838)


Description	HPV18-E6/E6 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from HPV18 E6 (75-158/158 aa)
Target	HPV18 E6 protein
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	11 kDa
Uniprot ID	P06463
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of paraffin...



IHC-P of human cervical carcinoma tissue...



Immunohistochemical analysis of formalin...