

Product Datasheet

Sphingosine 1 Phosphate Receptor 1 CT peptide (orb108394)

Description	Sphingosine 1 Phosphate Receptor 1 (EDG-1) CT blocking peptide
Conjugation	Unconjugated
Tested Applications	FA
Preservatives	Provided in solution at 1 mg/ml in phosphate buffered saline.
Form/Appearance	Provided in solution at 1 mg/ml in phosphate buffered saline.
Concentration	See vial for concentration
Storage	Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles
Note	For research use only
Application notes	For use with EDG1 (S1P1) (CT) polyclonal antibody (Cat. No. X1093P).
Hazard Information	This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.
Expiration Date	12 months from date of receipt.