

Product Datasheet

SARS-CoV-2 Spike Protein S1 antibody (orb1084011)

Description

SARS-CoV-2 Spike Protein S1 antibody

Species/Host

Rabbit

Conjugation

Unconjugated

Tested Applications

ELISA

Immunogen

KLH conjugated synthetic peptide derived from 2019-nCoV Spike Protein S (21-120/1273 aa)

Target

SARS-CoV-2 Spike Protein S1

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Liquid

Concentration

1mg/ml

Storage

Store at -20°C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4°C.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

140/65 kDa

Uniprot ID
P0DTC2
Expiration Date

12 months from date of receipt.