

## Product Datasheet

# SARS-CoV-2 (2019-nCoV) Nucleocapsid antibody (orb1084018)

**Description** SARS-CoV-2 (2019-nCoV) Nucleocapsid antibody

**Species/Host** Mouse

**Reactivity** Virus

**Conjugation** Unconjugated

**Tested Applications** WB

**Immunogen** Recombinant SARS-CoV-2 (2019-nCoV) Nucleocapsid protein (full length)

**Target** SARS-CoV-2 (2019-nCoV) Nucleocapsid

**Preservatives** 0.01M PBS(pH7.4) with 0.03% Proclin300

**Form/Appearance** Liquid

**Concentration** ≥1mg/ml

**Storage** Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

**Note** For research use only


**Isotype** IgG1

**Clonality** Monoclonal

**Uniprot ID** [P0DTC2](#)

**Dilution Range** WB=1:500-2000

**Expiration Date** 12 months from date of receipt.



Western blot analysis using Anti-SARS-Co...