

Product Datasheet

Sphingosine 1-Phosphate Phosphatase 1 peptide (orb108402)

Description	Sphingosine 1-Phosphate Phosphatase 1 blocking peptide
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FA
Preservatives	Provided as solution in phosphate buffered saline with 0.08% sodium azide
Form/Appearance	Provided as solution in phosphate buffered saline with 0.08% sodium azide
Concentration	See vial for concentration
Storage	Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles
Note	For research use only
Application notes	For use with Sphingosine 1-Phosphate Phosphatase 1 polyclonal antibody (Cat. No. X1636P).
Isotype	IgG
Uniprot ID	Q9BX95
Hazard Information	This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.
Expiration Date	12 months from date of receipt.