

Product Datasheet

human CD11b antibody (orb1084046)

Description	human CD11b antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	Recombinant human CD11b
Target	human CD11b
Preservatives	0.01M TBS(pH=7.4)
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1, k
Clonality	Monoclonal
MW	125 kDa
Uniprot ID	P11215
Dilution Range	Flow-Cyt=1ug/Test
Expiration Date	12 months from date of receipt.